

BST Thick-Film Sensors

Art. 100 000000

2 Electrodes

Substrate: Al₂O₃ ceramic

Dimensions (L x W x T): 25.4 x 7.0 x 0.635 mm

WE – Carbon (diameter 1.0 mm)

REF – Ag/AgCl

Isolator paste: low temperature



Art. 154 000000

2 Electrodes

Substrate: Al₂O₃ ceramic

Dimensions (L x W x T): 25.4 x 7.0 x 0.635 mm

WE – Gold (diameter 1.0 mm)

REF – Ag + Ag/AgCl

Isolator paste: high temperature



Art. 120 000000

3 Electrodes

Substrate: Al₂O₃ ceramic

Dimensions (L x W x T): 25.4 x 7.0 x 0.635 mm

WE – Platinum (diameter 1.0 mm)

CE – Platinum

REF – Ag

Isolator paste: high temperature



Art. 1212 000000

3 Electrodes

Substrate: Al₂O₃ ceramic

Dimensions (L x W x T): 25.4 x 7.0 x 0.635 mm

WE – Gold (diameter 1.0 mm)

CE – Platinum

REF – Ag + Ag/AgCl

Isolator paste: high temperature

